Build-It™ FR

Core Material

English

INSTRUCTIONS FOR USE

Français

MODE D'EMPLOI

French

Deutsch

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German

Italiano

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Italian

Español

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Spanish

Čeština

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Czech

English

INSTRUCTIONS FOR USE

Build-It™ F.R.™ Core Material contains all of the features of an ideal core build-up material. Formulated to flow and stack without slumping, the material's unique viscosity makes it the perfect core material as well as the ideal post cement. Specially treated glass fibers give Build-It F.R. Core Material a compressive strength of 329 MPa for a strong and reliable foundation for your indirect restorations. Designed to offer optimal versatility and performance, Build-It F.R. Core Material self-cures in minutes or light-cures in seconds to a rock-hard consistency that cuts like dentin without ditching. Choose from three esthetic dual-cure shades (A2, A3, & Opaqueous White) or two self-cure contrast shades (Gold, Blue) in either 25ml auto-mix cartridges or 4ml Mini Mix™ syringes.

Indications: Build-It FR is a fiber reinforced, dual cure post and core build-up resin composite material used for restoring a crown or missing tooth structure.

LOADING THE CARTRIDGE DISPENSER GUN

1. Push up the small lever in the back of the Dispenser Gun and make sure the slide (serrated side facing down) is inserted completely into the Dispenser Gun.
2. Lift the plastic hinged cover on the top of the gun and insert the cartridge. NOTE: A notch on the bottom of the cartridge aligns with the notch on the Dispenser Gun. Push down firmly on the hinged cover to lock the cartridge in place.
3. Before attaching the mixing tip, remove the cartridge cap by turning it 90° counter clockwise, pull and remove. Dispense out a pea size amount of material and discard.
4. Line up the notch on the mixing tip with the corresponding notch on the cartridge flange. Push the mixing tip onto the cartridge and turn it clockwise 90° until it locks into position. If desired, an intra-oral tip can be inserted into the open end of the mixing tip.
5. Pull the trigger of the Dispensing Gun with even pressure to dispense Build-It F.R. material.
6. **TO REMOVE THE CARTRIDGE:** Push up the lever on the back of the Dispenser Gun and pull the slide back all the way. Lift the hinged cover up on the top of the gun and remove the cartridge.

7. Leave the mixing tip on the cartridge until next use. The tip serves to seal the cartridge.

USING THE MINI MIX™ SYRINGE

1. Remove the syringe cap by turning it 90° counter clockwise, pull and remove. Express core material prior to placing the mixing tip onto the syringe to ensure that both sides flow evenly.
2. Line up the notch on the mixing tip with the corresponding notch on the cartridge flange. Push the mixing tip onto the cartridge and turn it clockwise 90° until it locks into position. If desired, an intra-oral tip can be inserted into the open end of the mixing tip.

LOADING THE AUTO-MIX SYRINGE DISPENSER GUN

1. Remove the plunger from the Auto-Mix syringe.
2. Push up the small lever in the back of the Auto-Mix Syringe Dispenser Gun and pull the plunger (serrated side facing down) all the way back.
3. Slide the Auto-Mix syringe into the dispenser gun until it locks into place.

IMPORTANT NOTICE:

If you are using a resin based temporization material or a resin based temporary cement (e.g.-TempSpan™ Temporary Crown & Bridge Material, TempSpan Temporary Cement) over the Build-It F.R. core, you must put a separating medium (petroleum or KY Jelly) on the core to prevent the temporary material from adhering to it.

CORE BUILD-UP

1. Isolate the tooth with a rubber dam.
2. Apply 37% Phosphoric Acid Etching Gel* to enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etch into eyes, onto soft tissues or skin.
3. Dry with clean, oil-free, air without desiccating the dentin. Do not allow etched surface to become contaminated with saliva or blood during the bonding procedure. NOTE: Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.
4. Rinse the tooth with copious amounts of water to remove the etchant. Blot or lightly blow dry for 2 seconds to remove excess moisture. Do not over-dry, leave the surface moist.
5. Apply 2 consecutive coats of Bond-1™

Primer/Adhesive within 10 seconds to the dentin and enamel. Gently air dry to remove the solvent. Do not over-dry; the surface should be shiny with resin.

6. Light cure for 10 seconds.
7. **A. If using a matrix band** - Place a matrix band around the prepared tooth. Place the end of the intraoral tip directly onto the surface of the preparation and dispense Build-It F.R. material. Build-It F.R. material has a minimum working time of 45 seconds and a set time within 4 minutes. Build-It F.R. self-cures within 4 minutes. Alternatively, you can light cure the lingual, facial and occlusal aspects of the A2, A3, or Opaqueous White shades for 40 seconds each. Blue and Gold are self-cure only and can be prepped within 4 minutes from the start of mix.
- B. If using a core form** - Place the mixing tip directly into the core form and dispense Build-It F.R. material. Build-It F.R. material has a minimum working time of 45 seconds and a set time within 4 minutes. Place the loaded core form onto the prepared tooth and hold in place until Build-It F.R. material is fully set. Alternatively, you can light cure the lingual, facial and occlusal aspects of the A2, A3, or Opaqueous White shades for 40 seconds each. Blue and Gold are self-cure only and can be prepped within 4 minutes from the start of mix.
8. Prepare and finish the core in the usual manner using diamond and carbide burs.

NOTE:

Light Curing A2, A3 and opaqueous white; Cure the buccal, lingual and occlusal surfaces for 40 seconds each.

POST CEMENTATION USING Build-It F.R.

1. Apply 37% Phosphoric Acid Etching Gel* to enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etch into eyes, onto soft tissues or skin.
2. Dry with clean, oil-free, air without desiccating the dentin. Do not allow etched surface to become contaminated with saliva or blood during the bonding procedure. NOTE: Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.
3. Flush the post hole with a copious amount of water and lightly blow dry. Use paper points to remove any excess moisture. Do not over-dry,

leave the surface moist.

- Dispense 1 drop of Bond-1 Primer/Adhesive via needle tip syringe into the post hole preparation. Allow it to sit for 20 seconds.
- Wick out any excess with an absorbent paper point. Do not leave any pools of bonding agent. Do not over-dry.
- Coat a Lentulo Spiral or paste filler with Build-It F.R. Core Material. Place the spiral into the post space and spin the Build-It F.R. Core Material along the canal walls while removing the spiral.
- Coat a FibreKor™ or FibreKleer™ Post with Build-It F.R. Core Build Material and insert it into the canal with gentle finger pressure, allowing excess core material to vent.
- Insert the post into the post space allowing the excess Build-It F.R. material to vent. Remove the excess Build-It F.R. material.
- Build-It F.R. material sets within 4 minutes in the post space. Allow the Build-It F.R. material to set before trimming the post. Alternatively, you can light cure the occlusal aspect for 40 seconds and then trim the post.
- Fabricate your core following the instructions above.

NOTE: Tip and applicators brushes are for single patient use only. This is to prevent cross-contamination between patients.

NOTE: Refer to outer package for expiration date. Do not use after expiration date.

***CAUTION:** 37% Phosphoric Acid Etching Gel contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water; get medical attention for eyes. Do not take internally.

CAUTION: Bond-1 Primer/ Adhesive contains acetone and is flammable and volatile. Keep away from sources of ignition and excessive heat.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Core Material

Rx Only

Refrigeration required.
Store at 2-8°C. (36-46°F)

LIMITED WARRANTY —

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

PENTRON

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / www.pentron.com

EC REP European Union Representative

SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111